visomat[®] comfort III

Gebrauchsanleitung (D2 - D27)

Instruction manual (GB28 - GB53)

Mode d'emploi (F54 - F79)

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CE 0123

Inhaltsangabe

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Contents

A	Product information/notes on use	30
	Device description/indicators	32
B	Abridged instructions on measurement	34
C	Operation of the device	
	1. Setting up the device	36
	2. Fitting the cuff	36
	3. Inflating the cuff and measuring blood pressure	38
	4. Irregular pulse waves	39
	5. Calling up the memory	40
	6. Deleting the memory	40
D	Important information	
	1. Basic rules for blood pressure measurement	41
	2. Other points to note	41
	Arm circumference	41
	Cardiac arrhythmia	41
	Frequency of pacemakers	42
	Pregnancy	42
	Medicines	42
	Intravenous supply systems	42
	Unusual results	42

E What you should know about blood pressure

- | | |
|---|----|
| 1. Systolic and diastolic blood pressures | 43 |
| 2. Why you have to measure different values | 43 |
| 3. Why measure regularly? | 44 |
| 4. What are normal blood pressure levels? | 44 |

F Maintenance

- | | |
|---------------------------|----|
| 1. Changing the batteries | 45 |
| 2. Care and maintenance | 45 |
| 3. Customer service | 46 |

G Appendix

- | | |
|---|----|
| What to do in the event of measuring errors and operating errors? | 46 |
| Specifications/accessories/spare parts | 48 |
| Explanation of symbols | 50 |
| Metrological inspection | 51 |
| Warranty | 52 |

These instructions are designed to help the user operate the digital visomat® comfort III blood pressure monitor safely and efficiently.

The monitor must be used according to the procedures contained in these instructions and for no other purpose

Read through the instructions carefully and familiarise yourself particularly with the section entitled «Basic rules for blood pressure measurement».

Notes on use

The visomat® comfort III is a digital automatic monitor for measuring blood pressure at the upper arm. It is based on the non-invasive analysis of the systolic and diastolic blood pressures at the upper arm and the measurement of adults' pulse rates, i.e. those of people 15 years of age and over. If the monitor is used on younger patients, the readings should be compared with doctors' measurements. Please call our Service phone line for more information.

Do not under any circumstances use the monitor on babies/infants.

The product is recommended to patients with unstable blood pressure or known hypertension, for home use and therapy support.

There may be a flaw in the accuracy of the measurement if absolute quiet is not observed while the measure-

Notes on use

ment is being taken or the conditions stated are outside those parameters laid down in the specifications. See technical data.

The cuff fits upper arms between 23 and 43 centimetres in circumference.

The device covers pressures up to 300 mmHg and a pulse rate of 40 to 160 beats per minute.


Method of measurement

This product uses the oscillometric method for measuring blood pressure and pulse rate. The cover is connected to the device and placed round the arm as shown on the diagram on the cuff. When the Start/Stop button is pressed, the device starts to build up pressure automatically.

Within a very short space of time the device records the small oscillations within the cuff, which occur owing to expansion and contraction of the arteries in the arm (pulsation). The automatic pressure buildup usually stops approx. 40 mmHg above the systolic value (fuzzy logic). Then the measuring process proper starts by letting air out of the cuff. The amplitude of each pressure wave is measured in millimetres of mercury (mmHg), converted and displayed on the LCD in the form of a digital value. Systole, diastole and pulse are established and after that the cuff is completely deflated.

The monitor offers an additional display for identifying irregular pulses. This can be an indication of an irregular heartbeat (arrhythmia) but also of unrest during the measurement (Irregular pulse waves).

A memory stores the last 30 results for the purposes of comparison.

Device description

A

- ❶ Start/Stop button
- ❷ Memory button
(memory recall button)
- ❸ Displays for systole, diastole,
pulse and indicators
- ❹ Air socket
- ❺ Socket for mains adaptor
- ❻ Battery compartment
- ❼ Cuff
- ❽ Mark for arm circumference
- ❾ Measurement arrow for arm circumference
- ❿ Artery marking

Indicators/symbols

Device measuring

Irregular pulse waves

Device inflating

Cuff deflating

Change batteries

*See special notes, P.45

Err

Measuring errors Err-1, -2, -3, 300

SYS

Systole

DIA

Diastole

PUL

Pulse rate

Memory symbol

Please read the detailed Instructions for Use before measuring for the first time (from page 36 onwards).

1. Insert batteries

2. Fit the cuff

Bare your left arm.



Place the cuff on the upper arm:

The bottom edge of the cuff should be 2-3 cm above the elbow and the hose should be at the middle of the inside of the elbow. Important!

The mark (approx. 4 cm bar) has to be directly above the artery running along the inside of the arm.

Close the cuff and connect it up to the device.

Stretch out your left arm and place it on a table, with the palm facing up. Keep your arm still and do not talk.

Abridged instructions on measurement

3. Inflate the cuff

Press the Start/Stop button. The device conducts a performance test and all the display elements light up.

The cuff is inflated to the necessary pressure.

The actual measuring process then begins. The „♥“ symbol on the display flashes to indicate the pulse rate and there is an audible beep. The measurement is taken, then the cuff is fully deflated. The end of measurement is indicated by a long beep.

The display shows the blood pressure levels of systole and diastole along with the pulse rate.

4. Switch off the device

Press the Start/Stop button. The figures measured are saved automatically.

1. Setting up the device

Place the enclosed batteries in the device if they have not already been inserted (see page 45, Changing batteries).


If you want to operate the device from the mains, insert the plug on the mains adaptor (see Accessories) into the socket on the right side of the device.

The batteries are switched off automatically.

2. Fitting the cuff


Bare your left upper arm. Slide the open cuff onto your upper arm until the bottom edge of the cuff is 2-3 cm above the inside of the elbow. When the cuff is used on the left arm, the air tube runs over the middle of the inside of the elbow to the monitor; in this way, the 4 cm long artery marking is located over the pulse position.

When the cuff is used on the right arm, turn the sleeve as far to the left as possible until the artery

Important:

Never use the cuff above a critical part of the body (e.g. wound, intravascular access).



Check whether the arrow on the cuff is within the white mark at the edge of the cuff (arm circumference 23-43 cm).

Operation of the device

marking is located on the pulse position. The air tube then runs along the inside of the upper arm.

The instructions for use apply by analogy. The cuff should fit in such a way that 2 fingers can pass between the arm and the cuff. For this purpose, bend your arm slightly before closing the Velcro fastener. The upper arm muscle is tensed slightly and the circumference of the arm increases slightly. Now pull the free end of the cuff tight and place it over the Velcro fastener firmly. Insert the cuff air connector into the air socket on the device as far as it will go. Lay your upper arm – stretched out and relaxed – on the table with the palm turned upwards. Keep the arm still and do not talk.

Important:

The cuff may only be inflated when it is attached to the arm.

Before measurement take a deep breath and relax. Do not move, do not talk and do not tense your arm during measurement. Breathe calmly.

3. Inflating the cuff and measuring

Press the Start/Stop button to switch the monitor on. The entire display lights up and this enables you to check that all the display elements are working correctly (fig. 1).

fig. 1

The automatically controlled inflation process starts after the monitor has calibrated itself against the ambient air pressure.

Then the measuring process proper starts with deflation. The „♥“ symbol on the display flashes and there is a beep for pulse rate (fig. 2).

fig. 2

The determined values for systole, diastole and pulse are shown in the display and at the same time the cuff is fully deflated.

A long beep signifies the end of the measuring process.

Note:

The inflation and measuring processes can be interrupted at any time by pressing the Start/Stop button.

Operation of the device

4. Irregular pulse waves

If the irregular pulse wave symbol flashes after a measurement, the monitor has registered an irregular pulse during the measurement. This can be set off by an irregular heartbeat (arrhythmia), disturbances caused by movement, talking or even breathing in deeply. Essentially, it involves the detection of an event which can unfavourably influence the measurement results. Measurement results accompanied by flashing pulse displays should therefore be viewed as critical and repeated under more favourable conditions.

Constant repetition of measurements can bring about congestion and thus lead to incorrect results. Between measurements take a break of at least 3-5 minutes.

5. Calling up the memory

Measured results are automatically stored in the memory. The memory can store up to 30 results and the average value.

When more than 30 measured values have been stored, the oldest value (No. 30) is deleted so that the most recent value (No. 1) can be recorded.

Press the memory button to recall data.

The average value of the stored results saved is displayed with "A" followed by the most recently measured value. If no value is stored, nothing is displayed.

Repeatedly press the button to call up the other measured values. The stored data are displayed for about 30 seconds. The monitor will then switch itself off.

6. Deleting the memory

To delete individual measurements, press the Memory button a number of times to call up the measurement required, press again and keep the button pressed for 8 - 10 seconds until the measurement has disappeared. To delete the entire memory, call up the mean, press the Memory button again and keep it pressed until the mean has disappeared (8-10 sec).

Press the Start/Stop button and switch the monitor off.

If the display flashes, the deletion process has not been completed and must be repeated.

Record the figures saved in your blood pressure pass. Over a lengthy period this will provide you with a better picture of your blood pressure situation, which you may wish to discuss with your doctor.

Important information

1. Basic rules for blood pressure measurement

- Measurement must be performed when you are at rest. Blood pressure is increased by walking and climbing stairs. Wait at least 5-10 minutes until your circulation has calmed down.
- Check your blood pressure twice a day – in the morning after getting up, and in the evening when you have relaxed after work.
- Blood pressure must be measured when you are seated. For measurement, sit down at your dining table (never at a low coffee table) in as relaxed and upright a posture as possible, support the arm and back, and do not cross your legs.
- The most reliable way of obtaining inaccurate figures is to perform measurement on the bare arm. Wearing light clothing, e.g. shirt or blouse, you can also measure through the sleeve.

If you are wearing tight-fitting clothes, do not constrict the arm by rolling up your sleeve but

simply remove the article of clothing.

- During measurement the arm and the cuff must remain still in order to prevent interfering pulses. Pressing the cuff and kinking the tube can affect the results.
- Always wait roughly 3 to 5 minutes between successive measurements to allow the artery to regenerate.

2. Other points to note

- Arm circumference
The cuff is designed for an upper arm circumference of 23-43 cm. A smaller cuff is available as an option for arm circumferences of 17 to 25 cm, see Accessories.
- Cardiac arrhythmia / Irregular pulse waves
If you suffer from cardiac arrhythmia. Therefore consult your doctor before measuring your blood

pressure. The monitor is also equipped with a pulse unrest detection facility. This indicates to users the irregular pulse waves signals during a measurement. These can be caused by an irregular heartbeat or unrest and disruptive movements during the measurement. In each case, clarify all possible instances of an irregular heartbeat with your doctor.

- Frequency of pacemakers

The display of pulse rate is not suitable for monitoring the frequency of pacemakers.

- Pregnancy

Blood pressure levels can be different as a result of pregnancy. Only perform blood pressure self-monitoring after consulting your doctor.

- Medicines

Never use this device without consulting your doctor first if you are undergoing dialysis treatment or taking anti-coagulant, anti-platelet medicines

or steroids. The use of a blood pressure monitor under such conditions can cause internal bleeding (possibly causing haematomas).

- Intravenous supply systems

Medical supply systems (e.g. infusions) may temporarily discontinue supply if the cuff is constricted.

- Unusual results

Any measurement is largely dependent on the circumstances during the measuring procedure. Before measurement people often fail to take time to calm their circulation adequately, use incorrect cuff sizes or they talk or move about during measurement. See also the table of measuring errors and operating errors on page 47.

Check the situation yourself if necessary. Have you met all the conditions? Take time for measurement – you are monitoring your health!

What you should know about blood pressure

If you would like to have specific information about your blood pressure, ask your doctor. Never draw any conclusions of your own on the basis of your results.

Self-measurement does not constitute treatment. Never change of your own accord the dose of medicines prescribed by your doctor.

The use of this device near mobile telephones, micro wave appliances and any other equipment producing strong electromagnetic fields can cause malfunctions.

Store the device outside children's reach. Never place the air hose round your neck because there is a risk of an accident due to strangulation.

1. Systolic and diastolic blood pressures

The cardiovascular system has the important function of supplying all the organs and tissues in the body with sufficient blood and removing metabolic products. The heart contracts and relaxes at a regular rate of about 60-80 times a minute. The pressure of the flowing blood on the arterial walls, which arises when the heart contracts, is termed systolic. The pressure in the ensuing relaxation phase when the heart is refilled with blood is termed diastolic. When you perform daily monitoring you measure both figures.

2. Why you have to measure different values

Our blood pressure responds to external and internal influences like a sensitive measuring instrument. It can be affected by even slight changes. That makes it easier to understand why figures measured at the doctor's or pharmacist's are often higher than those you

What you should know about blood pressure

measure at home in your accustomed environment. Changes in weather, changes in climate, physical work and emotional stress can also have an effect.

3. Why measure blood pressure regularly?

The time of day also has an influence on the level of blood pressure. During the day the figures are usually higher than during periods of rest at night. For this reason, one-off measurements and irregular measurements will tell you little about your actual blood pressure. Reliable assessment is only possible if individual measurements are taken regularly.

Discuss the results with your doctor.

4. What are normal blood pressure levels

The World Health Organization (WHO) has compiled the following overview for assessing blood pressure levels:

Normal pressure	Threshold	High pressure
Systolic pressure <119	Systolic pressure 120–139	Systolic pressure >140* and /or**
Diastolic pressure < 79	Diastolic pressure 80–89	Diastolic pressure >90

* mmHg = millimetres of hydrargyrum (mercury)

** It is sufficient if one of the figures is too high

WHO 2003

Maintenance

1. Changing batteries

- Open the battery compartment:
Remove the cover of the battery compartment on the underside of the device.
- Insert batteries:
Take the old batteries out of the device and insert the new ones. Make sure the polarity is correct (there is a mark inside the battery compartment).
- Close the battery compartment:
Close the battery compartment by replacing the cover.
- Please note:
Change the batteries:
 - if after the segment check the flat battery symbol “” appears in the display.
 - if the display fails to light up despite your having pressed the Start button.

Always change all 4 batteries.

Only use long-life alkaline batteries (LR 6). Carbon zinc batteries (R 6) have to be changed frequently. Do not use rechargeable batteries. If the device is going to be out of use for any lengthy period, please remove the batteries. This prevents leakage. Any battery can leak.

Batteries and appliances should not be disposed of with domestic waste but should be handed in at appropriate collection and disposal points.

Caution!

Do not use rechargeable batteries.

2. Care and maintenance

The device contains sensitive components and should be protected against substantial

temperature fluctuations, humidity, dust and direct sunlight. The device is not impact-resistant or shockproof.

Protect against moisture. The device is not waterproof. Only use a soft, dry cloth to clean the device. Never use petrol, thinners or any other strong solvents. Do not scrub the cuff or wash it in the washing machine. Use a synthetic detergent and rub the surface gently.

Make sure no liquid gets into the air tube.

Allow to dry thoroughly.

If you have to disinfect the device, we recommend using Descosept AF, Dr. Schumacher GmbH, Melsungen. It can be used for fast wipe-down and spray disinfection.

The cuff must not be washed or dry-cleaned.

If the device has been opened, it must be subjected to a metrological inspection by an authorised institution.

3. Customer service

The device may only be repaired by the manufacturer or an authorised service centre.

Please contact:

UEBE Medical GmbH

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97877 Wertheim, Germany

Tel.: 09342/924040

Fax: 09342/924080

E-mail: info@uebe.com

Internet: www.uebe.com

What to do in the event of measuring errors or operating errors?

People with an irregular heartbeat, vascular constrictions, arteriosclerosis in their extremities, and diabetics should first consult their doctors before measuring their own blood pressures as in such cases different readings may be obtained.

Appendix

Malfunc- tions	Possible reasons	Corrective action
Err -1	Movement of the arm or body during the measurement. Do not talk during the measurement.	Repeat the measurement after 3-5 minutes. Don't move the arm and don't talk. In case of recurrence please contact our after-sales service.
Err -2	Pump problem	Check whether the air connection is correctly inserted. Is the cuff fitted too loosely? Repeat measurement. In case of recurrence please contact our after-sales service.
Err - 3	Deflation rate too high or too low, valve problem	In case of recurrence please contact our after-sales service.
Err - 300	Upper pressure limit of 300 mmHg has been exceeded.	Repeat the measurement after 3-5 minutes rest
	Irregular heartbeat, disruptive movements, shaking, wobbling (objects), breathing in deeply, etc.	Repeat measurement after 3-5 minutes rest. Possible factors (depending on the severity of the arrhythmia) affecting the measurement results should be discussed with your doctor.
Levels too high	Was the necessary rest observed before measurement?	Repeat measurement after a break of approx. 3-5 min. Do not move your arm and do not talk.
Unusual levels	Movement or talking during measurement. Resting time not observed, feet possibly crossed, refrain from smoking or drinking coffee at least 30 minutes before measurement.	Check conditions and repeat measurement.

Specifications/accessories/spare parts

- Type of device:
visomat® comfort III - automatic digital device with electric pump for blood pressure monitoring on the upper arm.
- The device bears a serial number SN which provides clear identification
- Device dimensions:
L= 140 mm x W = 130 mm x H = 70 mm
- Weight:
315 g not including batteries
- Measuring system:
Oscillometric determination of systole, diastole and pulse
- Reference procedure of clinical trial:
Auscultatory measurement
- Display:
LCD display (liquid crystal display)
- Memory:
30 measurements (saved automatically) and average value
- Pressure display range:
0–300 mmHg
- Measuring range:
Systolic: 50 - 250 mmHg
Diastolic: 40 - 150 mmHg
Pulse measurement: 40–160 Puls/min.
- Limits of error of static pressure display:
conforms to EN 1060 Part 3
Pressure display: ± 3 mmHg
Pulse measurement: $\pm 5\%$
- Cuff:
Cuff Type UPW (for arm circumference 23 – 43 cm)

Appendix

or also the optional
cuff Type US (for arm circumference 17-25 cm)

- Deflation valve:
Electronically controlled
- Inflation pressure:
Fuzzy logic controlled
- Deflation rate:
On average 4.5 mmHg/sec $\pm$ 1.5
- Batteries:
– battery type:
4 x 1.5 V, LR 6, size AA
alkaline or lithium (FR 6)
- Automatic shutdown:
approx. 3 minutes after the end of measurement
- Operating conditions:
Ambient temperature 10 – 40°C
Rel. humidity up to 85%

- Storage and transport conditions:
Ambient temperature -20 until +50°C
Rel. humidity up to 85%

Important:

The device may only be operated with genuine spare parts. If damage is caused by other makes of accessories, the warranty is invalidated!

Accessories

You can obtain the following genuine spare parts and accessories from specialist dealers:

- Cuff UPW (for arm circumference 23-43 cm)
Art. no. 2403001
- Cuff US (for arm circumference 17-25 cm)
Art. no. 2403005
- Mains adaptor A1
Art. no. 2401020

Subject to technical change without notice.

Explanation of symbols

This blood pressure monitor conforms to Council Directive 93/42/EEC of 14 June 1993 on Medical Devices and bears the mark CE 0123 (TÜV Süd Product Service). Devices with CE marks are quality-controlled in accordance with this Directive and provide a higher level of accuracy than previous calibration.

Degree of protection against electric shock:

TYP BF

Observe Instructions for Use

Batteries and appliances must not be disposed of with domestic waste but should be handed in at the appropriate collection and disposal points.

The requirements of the following standards are met:

- EN 1060-1: 1995
Non-invasive sphygmomanometers
Part 1: General Requirements
- EN 1060-3: 1997
Non-invasive sphygmomanometers
Part 3: Additional Requirements on electro-mechanical blood-pressure measuring devices
- EN 1060-4: 2004
Non-invasive sphygmomanometers
Part 4: Test procedure for establishing accuracy
- ESH (European Society of Hypertension)
The monitor satisfies the clinical testing requirements of the «International Protocol for validation of blood pressure measuring devices in adults, 2002».

Appendix

- EN 60601-1: 1990
 - + A1: 1993
 - + A2: 1995
 - 1. Medical electrical devices
 - 2. General commitment to safety
- IEC 60601-1-2: 2001 +A1:2004
 - Medical electric devices
 - Electromagnetic compatibility
- This product conforms to Council Directive 93/42/EEC of 14 June 1993 on Medical Devices.

Manufacturer:
UEBE Medical GmbH
Zum Ottersberg 9
97877 Wertheim
Germany

Important:

No modifications, e.g. opening the device (except to change the batteries), may be made to this device without the manufacturer's prior permission.

Metrological inspection

(previously calibration)

Generally, a metrological inspection is to be recommended at intervals of 2 years. Specialist and professional operators are referred to the relevant directives applicable in their individual countries. The inspection can be performed either by UEBE Medcial GmbH, an authority responsible for metrology or authorised maintenance service centres.

Competent authorities and authorised maintenance service centres can obtain a «Test Procedure for Metrological Inspections» from the manufacturer on request.

Warranty

The device was manufactured and tested with the utmost care. If it should nevertheless have defects, we give a warranty on the following terms:

1. During the term of the warranty of 3 years as of the date of purchase we agree to remedy such defects, at our discretion and at our expense, by repair (after the product has been returned to our factory) or in the form of a replacement by a non-defective device.
2. The warranty does not cover normal use of wearing parts or any damage which arises as a result of non-observance of the Instructions for Use, improper handling (e.g. unsuitable power sources, breakage, leaking batteries) and/or dismantling of the device by the purchaser. In addition, the warranty does not cover compensation claims filed with us.
3. Claims under the warranty can only be asserted within the warranty period and on submission of the original sales receipt. In the event of a warranty claim the device must be returned to the following address, together with the receipt and description of complaint:
UEBE Medical GmbH,
Service-Center,
Zum Schlag 18
97877 Wertheim, Germany.
4. The purchaser's contractual defect rights vis-a-vis the seller in accordance with § 437 of the German Civil Code are not restricted by the warranty.

Important:

In the event of a warranty claim please always include the receipt.

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© Copyright UEBE Medical GmbH